


ARS-IMAGO WB WASHING BATH




ARS-IMAGO WB - WASHING BATH is a concentrated liquid wetting agent to be used as the last bath after the final rinse. The treatment makes the film anti-static and prevents limestone traces on the surface. It also has a stabilising effect and reduces the bending of the films (in particular the ones with a PET component). It lowers the drying time.

Do not rinse the film at the end of the treatment.

PREPARATION

Dilute the concentrated liquid with running water before use, at the recommended dilution rate of 1+99.



**990ml WATER + 10 ml CHEMICAL
= 1L WORKING SOLUTION**

The final volume of the working solution depends on the size and kind of tank, the number of films/spirals and the size of trays that are being handled.

AGITATION

Agitate very gently without overturning the tank and remove the foam that may appear on the surface. To prevent the wetting agent from leaving traces on the spirals and contaminating subsequent treatments, it is advisable to immerse the film manually in the solution in a tray or jug.

TREATMENT TIME

Recommended treatment time: 30-60 seconds.

LIFE OF THE PRODUCT

The concentrated solution in its original packaging has no expiration date. It is recommended to use the diluted solution within 2 weeks after preparation.

KEEP OUT OF REACH OF CHILDREN

